


# **D-Stat® Dry Wrap**

Hemostatic Bandage

Topical Hemostat



***For the control of bleeding  
around indwelling lines***

- Pre-cut for easy placement around indwelling lines
- Contains the power of thrombin to facilitate hemostasis
- Clear bandage included for continuous site monitoring



# D-Stat® Dry Wrap

Hemostatic Bandage Topical Hemostat

Model	Description	Box Quantity
3015	D-Stat Dry Wrap Hemostatic Bandage	10*

\* Devices packaged in quantities of 10 units per box. Individual units are not sold separately.

The D-Stat Dry Wrap hemostatic bandage is applied topically as an adjunct to manual compression and is indicated for the control of surface bleeding from vascular access sites and percutaneous catheters or tubes.



*D-Stat Dry Wrap hemostatic bandage is part of Vascular Solutions' D-Stat line of hemostatic products, including D-Stat Flowable hemostat, D-Stat Dry hemostatic bandage, D-Stat Rad-Band topical hemostat, D-Stat 2 Dry hemostatic bandage and D-Stat Clamp Accessory.\*\**


The use of topical bovine thrombin preparations has occasionally been associated with abnormalities in hemostasis ranging from asymptomatic alterations in laboratory determinations, such as prothrombin time (PT) and partial thromboplastin time (PTT), to severe bleeding or thrombosis which rarely have been fatal. These hemostatic effects appear to be related to the formation of antibodies against bovine thrombin and/or factor V which in some cases may cross react with human factor V, potentially resulting in factor V deficiency. Repeated clinical applications of topical bovine thrombin increase the likelihood that antibodies against thrombin and/or factor V may be formed. Consultation with an expert in coagulation disorders is recommended if a patient exhibits abnormal coagulation laboratory values, abnormal bleeding, or abnormal thrombosis following the use of topical thrombin. Any interventions should consider the immunologic basis of this condition. Patients with antibodies to bovine thrombin preparations should not be re-exposed to these products.

The D-Stat Dry Wrap is contraindicated in persons with known sensitivity to bovine-derived materials. Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.  
CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

\*\* These D-Stat products are available in the United States. Not all products are available in all countries. Please contact the Vascular Solutions distributor or direct sales office for product availability.

U.S. and International patents pending.  
D-Stat is a registered trademark of Vascular Solutions, Inc.

©2009 Vascular Solutions, Inc. All rights reserved.  
ML2034 Rev. B 10/09



**Vascular  
SOLUTIONS**

**Vascular Solutions, Inc.**  
6464 Sycamore Court  
Minneapolis, Minnesota 55369 USA  
Customer Service: 888.240.6001  
[www.vascularsolutions.com](http://www.vascularsolutions.com)