

PRODUCT LAUNCH ANNOUNCEMENT



Monday, May 4, 2009

Contacts: Joan Will, Director of Marketing Communications, jwill@vascularsolutions.com
Vascular Solutions, Inc. 763.656.4300

Vascular Solutions Launches Amplatz SST™ guidewires

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. recently launched the Amplatz SST line of guidewires, a line of 0.035" and 0.038" guidewires designed for extra support during catheter placement in diagnostic angiography and interventional procedures.

Amplatz SST guidewires have an inner-core construction which provides added strength for increased support, a radiopaque tip for enhanced visibility and PTFE coating for enhanced device tracking. The Amplatz SST guidewire is available in straight and modified J configurations and in 0.035" and 0.038" diameters. Amplatz SST guidewires are currently available in the United States.



Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit www.vascularsolutions.com.